

ZFP36-zinc finger protein 36, C3H type, homolog (mouse) Gene

Reference: PTXBC009693

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ZFP36 ring finger protein; mRNA decay activator protein ZFP36; G0S24; GOS24; NUP475; RNF162A; TIS11; TTP; zfp-36; G0/G1 switch regulatory protein 24; growth factor-inducible nuclear protein NUP475; tristetrapoline; zinc finger protein 36 homolog; zinc fi		

For research use only.

GENE DESCRIPTION	zinc finger protein 36, C3H type, homolog (mouse)
GENE ID	7538
ORIGIN SPECIES	Human
ACCESSIONS	BC009693
VECTOR	pDONR223
NCBI SYMBOL	ZFP36
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatctgactgccatctacagagccctctgtcgtgagccctgacgtgccctgacatggaggagctgagtcagccaggtggggctctcgggacccctggagcctgagccctccgactccagccctgtgggtcactcccgctgctggccgtccacca
gcctagtgaggggccgagctgtggtgggtgccccacccctggcttcgacccgtggtcccccgctgggcccctgagctgtcaccctcaccacttgcgcaactccaccacccctcgcgtacaagactgagctatgtcggacttctcagagagtggcgctgccgc
tacggggccaagtgcagtttgccatggcctggcgagctgcccaggaatgcaccccaatacaagacgggaactctgtcacaagtctacctcagggcgctgcccctacggctctcgtgccacttcacacacccctagcgaagacctggcgcccgccacccctc
ctgtgcttcgagagcagcagcttctccgctgccccttgccgcccggacctcaccacccagggcctggcgcccttcctgtcctcagctccttctgcctccagctcccccaccacacccctggggaccttcacgtcaccctctgcttctctgctgccctggcaccctc
tggctcgaagagacccacccagctgtgttgcctctgcccgaaggccactctatcagcgtctgggggcccctgggtggcctgttgggacccctctgtacagctccctgggatccgaccctgatgaatatgccagcagcgagcagcctgggggctctgacttccgtcttcg
aggcgggagttttgcaccacccagccctggcagccccggcgactccccatcttcaatgcacatctgtttctgagtga