

G0S2-G0/G1switch 2 Gene

Reference: PTXBC009694

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	G0/G1 switch protein 2; G0/G1 switch regulatory protein 2; G0/G1switch 2; G0/G1 switch 2		

For research use only.

GENE DESCRIPTION	G0/G1switch 2
GENE ID	50486
ORIGIN SPECIES	Human
ACCESSIONS	BC009694
VECTOR	pDONR223
NCBI SYMBOL	G0S2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaacggtccaggagctgatccccctggccaaggagatgatggccagaagcgcaaggggaagatggtgaagctgtacgtgctgggcagcgtgctggccctcttcggcgtggtgctcgccctgatggagactgtgtgcagcccttcacggcccgacgtctcgggacca
ggagcgagccgtggcgagctgaggccgacctggagcgacaggctctccagaagcaagccctgcaggagaaagcaagcagcaggacacggctcgcggcgccctgtccaaccggcagcagcctcctag