

IFITM2-interferon induced transmembrane protein 2 (1-8D) Gene

Reference: PTXBC009696

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	1-8D; DSPA2c; interferon-induced transmembrane protein 2; dispanin subfamily A member 2c; interferon-inducible protein 1-8D; interferon induced transmembrane protein 2		

For research use only.

GENE DESCRIPTION	interferon induced transmembrane protein 2 (1-8D)
GENE ID	10581
ORIGIN SPECIES	Human
ACCESSIONS	BC009696
VECTOR	pDONR223
NCBI SYMBOL	IFITM2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaccacattgtgcaaaccttctctctgtcaacagcggccagcctcccaactacgagatgctcaaggaggagcagggaagtggctatgctggggcgccccacaacctgctccccgacgtccaccgtgatccacatccgcagcgagacacctgctgacctgacctgtgtcc
ctgttcaacacctcttcatgaacacctgctgctgggtctcatgcatctcctactccgtgaagtctaggacaggaagatggttggcgacgtgaccggggccagccctatgctccaccgccaaagtgcctgaacatctgggacctgatitgggcatcttcatgacattctgctgcat
catcccagtggtgtgtccaggccagcgatag