

APOC1-apolipoprotein C-I Gene

Reference: PTXBC009698

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Apo-CI; ApoC-I; apo-CIB; apoC-IB; apolipoprotein C-I; apolipoprotein C1		
For research use only.			

GENE DESCRIPTION	apolipoprotein C-I
GENE ID	341
ORIGIN SPECIES	Human
ACCESSIONS	BC009698
VECTOR	pDONR223
NCBI SYMBOL	APOC1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggctctctctgctccccgtcctgggtgtgtgtctgtcgtcgtcttgaaggccagccccagccaggggacccagacgtctccagtccttggataagctgaaggagtttgaacacactggaggacaaggctcgggaactcatcagccgcatcaaacagagtgaacttt
ctgccaagatcgggagtggttttcagagacatttcagaaagtgaaggagaaactcaagattgactcatga