

MT1E-metallothionein 1E Gene

Reference: PTXBC009699

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MT-1E; MT-IE; MT1; MTD; metallothionein-1E; metallothionein 1E (functional); metallothionein D; metallothionein-IE; metallothionein 1E		

For research use only.

GENE DESCRIPTION	metallothionein 1E
GENE ID	4493
ORIGIN SPECIES	Human
ACCESSIONS	BC009699
VECTOR	pDONR223
NCBI SYMBOL	MT1E
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacccaactgctcttgcgcactgggtgctctgcacgtgcgcggctcctgcaagtgcaagagtgcaaatgcacctctgcaagaagagctgctgttctgctgccccgtgggctgtgccaagtgtgccagggtgcgtctgcaaggggcatcggagaagtcagctgctgtgcctga