



Reference: PTXBC009701

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PAD-H19; PAD2; PDI2; protein-arginine deiminase type-2; peptidyl arginine deiminase, type II; protein-arginine deiminase type II; peptidyl arginine deiminase 2		

For research use only.

GENE DESCRIPTION	peptidyl arginine deiminase, type II
GENE ID	11240
ORIGIN SPECIES	Human
ACCESSIONS	BC009701
VECTOR	pDONR223
NCBI SYMBOL	PADI2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1