

GNG11-guanine nucleotide binding protein (G protein), gamma 11 Gene

Reference: PTXBC009709

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GNGT11; guanine nucleotide-binding protein G(I)/G(S)/G(O) subunit gamma-11; G protein gamma-11 subunit; guanine nucleotide binding protein (G protein), gamma 11; guanine nucleotide-binding protein G(I)/G(S)/G(O) gamma-11 subunit; G protein subunit gamma 1		

For research use only.

GENE DESCRIPTION	guanine nucleotide binding protein (G protein), gamma 11
GENE ID	2791
ORIGIN SPECIES	Human
ACCESSIONS	BC009709
VECTOR	pDONR223
NCBI SYMBOL	GNG11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcctgcccttcacatcgaagattgccagagaaggaagaaactgaaatggaagtgagcagcttcgcaagaagtggaagtcagagacaacaagtgctaaatgttctgaagaaataaagaactatattgaagaacgttctggagaggatcctctagtaaagggaattccagaagacaa
gaaccctttaagaaaaagcagctgtgtttattcataa