

GOSR2-golgi SNAP receptor complex member 2 Gene

Reference: PTXBC009710

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Bos1; EPM6; GS27; Golgi SNAP receptor complex member 2; 27 kDa Golgi SNARE protein; membrin		

For research use only.

GENE DESCRIPTION	golgi SNAP receptor complex member 2
GENE ID	9570
ORIGIN SPECIES	Human
ACCESSIONS	BC009710
VECTOR	pDONR223
NCBI SYMBOL	GOSR2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatccctgttccagcaaacgcacagcagggtccagagatccagtcgtgcatgggacgcctggagacggcagacacagcagtcgtgcatagtagaaaacgaaatccaagcaagcatagaccagatattcagccgtctagaacgtctggagattttgcccagcaaggagcccccct
aacaaaaggcaaatgccagactcgggtgaccagttaaagtatgatgtccagcacctgcagactgcgctcagaaactccagcatcggcgccatgcaaggagcagcagagagacagcgagaagagcttctgtctgaacctccaccactaacgactctgacaccaccataccaat
ggacgaatcactgcagtttaactctccctccagaaagttcacaacggcatggatgacctcatctttagatgggcacaataatttagatggactgaggacccagagactgacctgaaggggactcagaagaagatccttgacattgccaacatgctggcctgtccaacacagtgatgcggct
catcgagaagcgggcttccaggaagactttatgataggtgggagctgctgacctgtgtggtcatgttctcgtgtgctgacgtacatga