

FKBP7-FK506 binding protein 7 Gene

Reference: PTXBC009711

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PPIase FKBP7; peptidyl-prolyl cis-trans isomerase FKBP7; FKBP23; PPIase; 23 kDa FK506-binding protein; 23 kDa FKBP; FK506-binding protein 23; FKBP-23; FKBP-7; rotamase; FK506 binding protein 7		

For research use only.

GENE DESCRIPTION	FK506 binding protein 7
GENE ID	51661
ORIGIN SPECIES	Human
ACCESSIONS	BC009711
VECTOR	pDONR223
NCBI SYMBOL	FKBP7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccaaaaccatgcatttcttattcagattcattgtttctttatctgtggggcccttttactgtctcagagacaaaagaaggagagcaccgaagaagtgaataagaagtttgcacgtccagaaaactgctctaagacaagcaagaaggagagactactaaatgccattatgacggct
acctggctaaagacggctcgaattctactgcagccggacacaaaatgaaggccaccccaatggtttgttcttgggtgttgggcaagtcataaaaggcctagacattgctatgacagatatgtgccctggagaaaagcgaagaagttagttatacccttcatttgcatacggaaagggaaggct
atgcagaaggcaagattccaccggatgctacattgttttgagattgaactttatgtgtgaccaaaggaccacggagcattgagacatttaacaaatagacatggacaatgacaggcagctctctaaagccgagataaacctctacttgcgaagggaatttgaaaaagatgagaagccac
gtgacaagtcatacagatgcagttttagaagatattttaagaagaatgacatgatggtgatggcttcttcccaaggatacaatgtatacacacacgatgaactatag