

ABCD3-ATP-binding cassette, sub-family D (ALD), member 3 Gene

Reference: PTXBC009712

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ABC43; CBAS5; PMP70; PXMP1; ZWS2; ATP-binding cassette sub-family D member 3; 70 kDa peroxisomal membrane protein; ATP-binding cassette, sub-family D (ALD), member 3; Peroxisomal membrane protein-1 (70kD); dJ824O18.1 (ATP-binding cassette, sub-family D (A		

For research use only.

GENE DESCRIPTION	ATP-binding cassette, sub-family D (ALD), member 3
GENE ID	5825
ORIGIN SPECIES	Human
ACCESSIONS	BC009712
VECTOR	pDONR223
NCBI SYMBOL	ABCD3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggccttcagcaagtacttgacggcgcaaaactcctcgtggctgggtccgcttctgctgctcctgctccacaagcggcgccgcgcctcggcctgcacggaagaaagtggaaaaccaccattacagaacaatgagaagaggggaaaaaggagcgagctgtgtg
gacaaggtgttttcaaggctcatacagattctgaaatcatgggtccctagaacattttgtaagagacaggttacttggtacttattgctgtatgctggtgtctgaacataatgtgatgtttggatgattcaaaatgggacactaatgaaagtgtatcattggtcgtagcagaaagatttcaa
gagatactactcaactcatcgtcctcatgccccttctctctggttaataactcttgaagtatgggttaaatgagcttaactgtgctccgagtaaggctcactaaatcctctatgaggagtatcttcaagctttcacatattataaaatggggaatctggacaacagaatagctaaccagacc
agctgcttacacaagatgtagaaaaatttgaacagtgtagtcgatctgtattcaatcttagtaagccatttttagacatagttttgtatatcttaagtaacgagtgaattggagctcaggtacttggaaaaatttggcattaa