

RAB39B-RAB39B, member RAS oncogene family Gene

Reference: PTXBC009714

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAB39B, member RAS oncogene family; MRX72; WSMN; ras-related protein Rab-39B		

For research use only.

GENE DESCRIPTION	RAB39B, member RAS oncogene family
GENE ID	116442
ORIGIN SPECIES	Human
ACCESSIONS	BC009714
VECTOR	pDONR223
NCBI SYMBOL	RAB39B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggccatctggctgtaccagttccgctcattgtcatcggggaattccacagtgggcaagtccctgacccgccgtccaccagggtcgtttgccagggtttctgacccaccgtgggggtgattttttctccgctggaggagatcgagcaggaaaacgcataagctccag
atctgggataccgcgggtcaagagaggttcagatccatcactcgcgcctactacaggaactcagtagtggtctctcttattgacattaccaaccgcagggtcctccagaatgtccatgagtggttagaagagaccaaagtacacgttcagccctaccaaatgtattgtctgtgggtcac
aagtgtgacctggatacacagaggcaagtgactcggccagggccgagaactggctgctgacacggcatgaagtacattgaaacgtcagcccagatgccattaatgtggagaaagccttcacagacctgacaagagacatatatgagctgggttaaaggggggagattacaatcca
ggagggtcgggaagggtgaagagtgattgtaccaaatgtgttctactcttcagaagaggtgtcaaatcagagaggagatgtttgtgctag