

CCL2-chemokine (C-C motif) ligand 2 Gene

Reference: PTXBC009716

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GDCF-2; HC11; HSMCR30; MCAF; MCP1; SCYA2; SMC-CF; C-C motif chemokine 2; chemokine (C-C motif) ligand 2; monocyte chemoattractant protein-1; monocyte chemotactic and activating factor; monocyte chemotactic protein 1; monocyte secretory protein JE; small i		

For research use only.

GENE DESCRIPTION	chemokine (C-C motif) ligand 2
GENE ID	6347
ORIGIN SPECIES	Human
ACCESSIONS	BC009716
VECTOR	pDONR223
NCBI SYMBOL	CCL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaagtctctgccgccctctgtgctgctcatagcagccaccttcattcccaagggtcgtcagccagatgcaatcaatgccccagtcacctgctgttataacttcaccaataggaagatctcagtcagaggctcgcgagctatagaagaatcaccagcagcaagtgtccaaa
gaagctgtgatcttcaagaccattgtggccaaggagatctgtgctgacccaagcagaagtgggtcaggattccatggaccacctggacaagcaaacccaaactccgaagactga