

ZMAT5-zinc finger, matrin type 5 Gene

Reference: PTXBC009717

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SNRNP20; U11/U12-20K; ZC3H19; zinc finger matrin-type protein 5; U11/U12 small nuclear ribonucleoprotein 20 kDa protein; U11/U12 snRNP 20 kDa protein; U11/U12 snRNP 20K; zinc finger CCCH-type containing 19; zinc finger matrin-type 5		

For research use only.

GENE DESCRIPTION	zinc finger, matrin type 5
GENE ID	55954
ORIGIN SPECIES	Human
ACCESSIONS	BC009717
VECTOR	pDONR223
NCBI SYMBOL	ZMAT5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggaagcgatacttctgtgactactgcgaccgtctctccaggacaacctccacaaccgcaagaagcactgaacgggctgcagcactcaaggccaagaaggtctgttacgacatgttcgagatgcagctgccatcttctggtgatgagcagaacaagcgccctgcaggaagtt
tctactgacaggccagtgcgactttggctccaactgcagattttccacatgtcagagcgagacctgcaggagctgagcatccaggtggaggaggagagggcgagccaggagtggtctactagatgctcctgagctccccgaggccatctggaggactggctggagaagagagccaa
gcggctgagctcagcccaagtagcagggtgaacccatcagaaccactgtctccagtaacccgtgggctggccaccagttcaggagctgcctccatccctgcgggcacccccacctgggggggtggcctctgcagccagagtcagttggggctga