

EDN1-endothelin 1 Gene

Reference: PTXBC009720

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ARCND3; ET1; HDLCQ7; PPET1; QME; endothelin-1; preproendothelin-1		
For research use only.			

GENE DESCRIPTION	endothelin 1
GENE ID	1906
ORIGIN SPECIES	Human
ACCESSIONS	BC009720
VECTOR	pDONR223
NCBI SYMBOL	EDN1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggattattgctcatgatttctctctgtgttggctgccaaggagctccagaacacagcagctcttaggcgctgagctcagcgcggtgggtgagaacggcgggagagaacccactccagtcaccctggcggctccggtccaagcgtgctcctgctcctgatggataaa
gagtgtgtctacttgcacacctggacatcatttgggtcaacactcccagcagcgtgttccgtatggacttggagccctaggtccaagagagccttggagaatttacttccacaaaggcaacagaccgtgagaatagatgccaatgtgctagccaaaagacaagaagtgtggaatttt
gccaagcaggaaaaaactcaggcgtgaagacattatggagaaagactggaataatcataagaaggaaaagactgttccaagcttgggaaaaagtgtattatcagcagtttagtgagaggaagaaaaatcagaagaagttcagagggaacacctaagacaaccaggtcggagacat
gagaaacagcgtcaaatcatctttcatgatccaagctgaaggaatccctccagagagcgttatgtgacccacaaccgagcacatttgtga