

IL1RN-interleukin 1 receptor antagonist Gene

Reference: PTXBC009745

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DIRA; ICIL-1RA; IL-1RN; IL-1ra; IL-1ra3; IL1F3; IL1RA; IRAP; MVCD4; interleukin-1 receptor antagonist protein; IL1 inhibitor; intracellular IL-1 receptor antagonist type II; intracellular interleukin-1 receptor antagonist (icIL-1ra); type II interleukin-1		

For research use only.

GENE DESCRIPTION	interleukin 1 receptor antagonist
GENE ID	3557
ORIGIN SPECIES	Human
ACCESSIONS	BC009745
VECTOR	pDONR223
NCBI SYMBOL	IL1RN
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctttagagacgatctccgacctctgggagaaaaatccagcaagatgcaagccttcagaatctgggatgttaaccagaagacctctctatctgaggaaacacactagttgccgatactgcaaggacaaatgtcaatttagaagaaaagatagatgtgtaccattgagcctcat
gctctgttctgggaatccatggagggaagatgtgctgtctgttcaagctctgtgatgagaccagactccagctggaggcagtttaacatcactgacctgagcgagaaacagaaagcaggacaagcgcttcgcttcacgctcagacagtggtcccccaccaccagttttgagtctccg
cctgccccggttggtcctctgcacagcgtggaagctgaccagccgctcagcctcaccaatagcctgacgaaggcgtcatggtcaccaaattctactccaggaggacgagtag