

MRPL48-mitochondrial ribosomal protein L48 Gene

Reference: PTXBC009764

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CGI-118; HSPC290; L48MT; MRP-L48; 39S ribosomal protein L48, mitochondrial; mitochondrial ribosomal protein L48		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein L48
GENE ID	51642
ORIGIN SPECIES	Human
ACCESSIONS	BC009764
VECTOR	pDONR223
NCBI SYMBOL	MRPL48
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagcggaaacctgggaaagtgctgtgcctgaggacaataccattttaagcaagcctttctcttaagggtttagaacttcaggagagaagcccatctattctgtaggtggaattctactaagtatcagtcggccctacaagacaagcccaccacggcattggaaagtacaagcatt
aattaaagcagaagagcccaagaagaaggaaggaaagtggaagtgagagccatttaattggggacagattatgaatatgggttttaaatattcatctgactgcatatgatatgaccctggcagagagttatgccagtagtttcacaaactctgcaactctctccattaaagtcaggaaa
gttatgcaatgccaaacaaacatagaagtggtgagtgacaggacaaaggcagcaaatgtctctggactcagtgcttaccacccatgagcagtggttcagatcagcgtttgagtgctacgtttgcagaaattttcttggaaataatccaaagcagttctctgaaggagtcagactgtc
agtgaaggagcacactgaagaagacttcaaggacgattcaagctcgaccagaactggaagaactgttggccaagttgaagtag