

C20orf30-chromosome 20 open reading frame 30 Gene

Reference: PTXBC009769

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0414 transmembrane protein C20orf30; C20orf30; HSPC274; dJ1116H23.2.1; transmembrane protein 230		

For research use only.

GENE DESCRIPTION	chromosome 20 open reading frame 30
GENE ID	29058
ORIGIN SPECIES	Human
ACCESSIONS	BC009769
VECTOR	pDONR223
NCBI SYMBOL	C20orf30
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatgccgtccctaccaacctggctactggaatcccagtagtaaaagtgaatattcaaggctctccagcacagacgatggctacattgacctcagtttaagaaacccctccaaagattccctataaggccatcgcaattgccactgtgctgtttttgattggcgctttctcattattatagg
ctccctcctgctgacaggtacatcagcaaaaggggggcagaccgggcccgtccagtgctgatcattggcattctggtgttctaccggattttaccacctgcgcattctactatgcaccaaaggctaccgtggttactcctatgatgacattccagacttggatgactag