

LYRM2-LYR motif containing 2 Gene

Reference: PTXBC009782

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DJ12208.2; LYR motif-containing protein 2; LYR motif containing 2		
For research use only.			

GENE DESCRIPTION	LYR motif containing 2
GENE ID	57226
ORIGIN SPECIES	Human
ACCESSIONS	BC009782
VECTOR	pDONR223
NCBI SYMBOL	LYRM2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgttcccgcttaccgccagcgacgctaacgttaagcagttcgtaagaaggcaacaagttctctctacagaaggattttgcaacaattcggcaagttccaaatgattctgatcgcaaatctgaaagattggcaagagaagaattcagaagaacaaaagtgccaccgaa
gaggatacaatccggatgatgattactcaaggcaatgacagctcaaggagttagaaaaaacacttgccttagcaaaatcttaa