

HTN3-histatin 3 Gene

Reference: PTXBC009791

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HIS2; HTN2; HTN5; histatin-3; basic histidine-rich protein; histatin-4; histatin-5; histatin-6; histatin-7; histatin-8; histatin-9; histidine-rich protein 3; hst; histatin 3		

For research use only.

GENE DESCRIPTION	histatin 3
GENE ID	3347
ORIGIN SPECIES	Human
ACCESSIONS	BC009791
VECTOR	pDONR223
NCBI SYMBOL	HTN3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagttttttgttttgccttaattcttgctctcatgcttccatgactggagctgattcacatgcaaagagacatcatgggtataaaagaaaattccatgaaaagcatcattcacatcgaggctatagatcaaattatctgtatgacaattga