

NDUFB6-NADH dehydrogenase (ubiquinone) 1 beta subcomplex, 6, 17kDa Gene

Reference: PTXBC009801

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	B17; NADH dehydrogenase [ubiquinone] 1 beta subcomplex subunit 6; CI-B17; NADH dehydrogenase (ubiquinone) 1 beta subcomplex, 6, 17kDa; NADH-ubiquinone oxidoreductase B17 subunit; NADH-ubiquinone oxidoreductase beta subunit, 6; complex I, mitochondrial res		

For research use only.

GENE DESCRIPTION	NADH dehydrogenase (ubiquinone) 1 beta subcomplex, 6, 17kDa
GENE ID	4712
ORIGIN SPECIES	Human
ACCESSIONS	BC009801
VECTOR	pDONR223
NCBI SYMBOL	NDUFB6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacgggtacactccgatgagaaactcggctgcagcagctcgagagctgagaaggcgatggctgaaggaccaggagctgagccctcgggagccggtgctgccccacagaagatgggacctatggagaaattctggaataaattttggagaataaattcccttgaggaaa
atggtccatgggtatatacaaaaagagtatcttttttcaatcatgtactgtacctgtctggattattcattattacatgaagtatcatgtttctgaaaaacatatggcatagttgaaaagaagtcagaataattccctggatgataaattctggagactggagaagtaattccacaaatgaagaattt
cctgatcaacatcaataa