

C11orf63-chromosome 11 open reading frame 63 Gene

Reference: PTXBC009820

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C11orf63; chromosome 11 open reading frame 63		
For research use only.			

GENE DESCRIPTION	chromosome 11 open reading frame 63
GENE ID	79864
ORIGIN SPECIES	Human
ACCESSIONS	BC009820
VECTOR	pDONR223
NCBI SYMBOL	C11orf63
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgccattctgagtttgatgatcgaatcggggcaacggtatggagcccgacgcttagacgaggaggaaagccctcgatggggaagcctgcacgagatggaagaggaaagcaagtggaagagcagctcagatggctcgcgagcaaaaccacatacctgggaccagggcgcca
ataacaggcaacaacaaatagaagacaataatcagacctccgctatgaccggaactggaagagtaagaaaggaggaaggcagctgctgctggaagcgttgccggagtcacggacagctcttagaaaaatctgcttggctccctctaccctccagagacgtcaatgggaact
ctccgggggaaaaggcgagcagaagaagagtcacagagtgacgctcttacttggtgagtgaaattttaagcccaactatgagcatggtgccgctgcagcaagccgttttcagagctgagcgacagtgacctggaggagaagtcgagcagccttctccgtacgtgaaagagctcaagt
tcacataacgaggttttctgcccgggacgacgtggccctggcgaaaggaaatccaacaacatttgggaaaaaacaagctcacttgggattacccacccgaaaacggagcttattctcaactcacaataaaaaaaggagggaatctcatccagaacaggtcattttgagagtaaat
ctctccacgtgatggattcaaaccttttggagccagtaa