

NUP160-nucleoporin 160kDa Gene

Reference: PTXBC009822

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	nucleoporin Nup160; nuclear pore complex protein Nup160; 160 kDa nucleoporin; nucleoporin 160kD; nucleoporin 160kDa; nucleoporin 160		

For research use only.

GENE DESCRIPTION	nucleoporin 160kDa
GENE ID	23279
ORIGIN SPECIES	Human
ACCESSIONS	BC009822
VECTOR	pDONR223
NCBI SYMBOL	NUP160
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcgggggagccctggaacggagcttcgtggagctaaagcggagctgagcgcgaaaggccgagggcacttcgggaattcacagtctgcagcattgggactgcaaatgccgtggcggcgcgtaaaatacagtgaaagcgcgggagcgtttactacgtggagagtggtcaa
gtgttctccgttaaccagaacagggttcattcattggaagacctctggagatacattggagctgatggaggagtcactggacataaatctgttgaataatgccattcgccataaattccaaatgagtggtttacgtggagggtttatgtctgagactcagaatcgtgtgataatcttgatgta
accaatcaaacagtgacaggttactttaccacaccctccggatgtaggagtgtaagttggcgaatgcaataatctttatftccaaattactctgggtgtcacaaatgtagtctggagcgtgtctttggaattgaagaaattggattctcgttatccctaccaagcactttgatag
ctaccgcttaaatga