



Reference: PTXBC009828

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LASS4; Trh1; ceramide synthase 4; LAG1 homolog, ceramide synthase 4; LAG1 longevity assurance homolog 4		

For research use only.

GENE DESCRIPTION	LAG1 homolog, ceramide synthase 4
GENE ID	79603
ORIGIN SPECIES	Human
ACCESSIONS	BC009828
VECTOR	pDONR223
NCBI SYMBOL	LASS4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]