

RAB25-RAB25, member RAS oncogene family Gene

Reference: PTXBC009831

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAB25, member RAS oncogene family; CATX-8; RAB11C; ras-related protein Rab-25		
For research use only.			

GENE DESCRIPTION	RAB25, member RAS oncogene family
GENE ID	57111
ORIGIN SPECIES	Human
ACCESSIONS	BC009831
VECTOR	pDONR223
NCBI SYMBOL	RAB25
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggaatggaactgaggaagattataactttgtctcaagtggtgctgacggcgaatcagggtgtggggaagaccaatctactctccgattcacgcgcaatgagttcagccacgacagccgaccaccatcggggttgagttctccaccgcactgtgatgtggcaccgctgctgt
caaggctcagatctgggacacagctggcctggagcgggtaccgagccatcacctggcgtactatcgtggtgagtgggggccctcctggtgtttgacctaaccaagcaccagacctatgctgtggtggagcgatggctgaaggagctctatgacctgctgaagccacgacgtcgtcat
gctcgtgggtaacaaaagtgacctcagccagccggggaagtggccactgaggaaggcccgatgttcgtgaaaacaatggactgctctcctggagacctcagccctggactctaccaatgttgagctagcctttgagactgtcctgaaagaatctttgcgaagggtgtccaagcagaga
cagaacagcatccggaccaatgccatcactctgggcagtgcccaggctggacaggagcctggccctggggaagaaggggcctgttgcacgcctctga