

SKP1-S-phase kinase-associated protein 1 Gene

Reference: PTXBC009839

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EMC19; OCP-II; OCP2; SKP1A; TCEB1L; p19A; S-phase kinase-associated protein 1; OCP-2; RNA polymerase II elongation factor-like protein OCP2; SIII; cyclin A/CDK2-associated p19; cyclin-A/CDK2-associated protein p19; organ of Corti protein 2; organ of Corti		

For research use only.

GENE DESCRIPTION	S-phase kinase-associated protein 1
GENE ID	6500
ORIGIN SPECIES	Human
ACCESSIONS	BC009839
VECTOR	pDONR223
NCBI SYMBOL	SKP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcctcaattaagttgcagagttctgatggagagatatgtgaagttgatgtgaaattgccaaacaatctgtgactattaagaccatgttggaagattgggaatggatgatgaaggagatgatgaccagttcctctaccaaatgtgaatgcagcaatattaaaaagggtcattcagtggtgca
cccaccacaaggatgacctctctcctctgaagatgatgagaacaagaaaaagcgaacagatgatccctgtttgggaccaagaattctgaaagttgaccaaggaacacttttgaactcattctggctgcaaaactcttagacatcaaaagtttgccttgatgttacatgcaagactgttgcc
aatatgatcaaggggaaaactcctgaggagattcgcaagaccttcaatcaaaaatgactttactgaagaggagggaagccaggtacgcaagagaaaccagtgggtggaagagaagtga