

HPCAL1-hippocalcin-like 1 Gene

Reference: PTXBC009846

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BDR1; HLP2; VILIP-3; hippocalcin-like protein 1; calcium-binding protein BDR-1; visinin-like protein 3; hippocalcin like 1		

For research use only.

GENE DESCRIPTION	hippocalcin-like 1
GENE ID	3241
ORIGIN SPECIES	Human
ACCESSIONS	BC009846
VECTOR	pDONR223
NCBI SYMBOL	HPCAL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggcaaacagacagcaagctgcggcccgaggtgctgcaggacctgcgggagaacacggagttcaccgaccacagctgcaggagtggtacaagggcttctcaaggactgccccaccggccacctgaccgtggacgagttcaagaagatctacccaactttccctacgg
cgacgttccaagttccgagcagctctccgacacttcgacaccaacggcgacggcaccatcgactccgggagttcatattgcgtgagcgtgacctcgcgggcaagctggagcagaagctcaagtgggccttcagcatgtacacctggacggcaacggctacatcagccgc
agcgagatgctggagatcgtgcaggccatctacaagatgggtcgtctgtgatgaagatgccggaggtgagtcaccccgagagaagcgcacagacaagatcttcaggcagatggacaccaacaatgacggcaaaactgtccttgaagaattcatcagaggtgccaagagcgaccct
ccatcgtccgctcgtgcagtcgaccccgacgtgccagtcagttctga