



Reference: PTXBC009894

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATPSK2; BCYM4; SK2; bifunctional 3'-phosphoadenosine 5'-phosphosulfate synthase 2; 3-prime-phosphoadenosine 5-prime-phosphosulfate synthase 2; ATP sulfurylase/APS kinase 2; ATP sulfurylase/adenosine 5'-phosphosulfate kinase; PAPS synthase 2; PAPS syntheta		

For research use only.

GENE DESCRIPTION	3'-phosphoadenosine 5'-phosphosulfate synthase 2
GENE ID	9060
ORIGIN SPECIES	Human
ACCESSIONS	BC009894
VECTOR	pDONR223
NCBI SYMBOL	PAPSS2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]