



Reference: PTXBC009920

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HNOEL-iso; OLF44; olfactomedin-like protein 3; olfactomedin like 3		
For research use only.			

GENE DESCRIPTION	olfactomedin-like 3
GENE ID	56944
ORIGIN SPECIES	Human
ACCESSIONS	BC009920
VECTOR	pDONR223
NCBI SYMBOL	OLFML3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]