

C12orf57-chromosome 12 open reading frame 57 Gene

Reference: PTXBC009925

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C10; GRCC10; protein C10; likely ortholog of mouse gene rich cluster, C10; chromosome 12 open reading frame 57		

For research use only.

GENE DESCRIPTION	chromosome 12 open reading frame 57
GENE ID	113246
ORIGIN SPECIES	Human
ACCESSIONS	BC009925
VECTOR	pDONR223
NCBI SYMBOL	C12orf57
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgtcgcctcgacccaacccggcgctgagcgtgagcaagcaagggtgctcgcggaggatgacaggcgttctccgccccggagaatgcagtcgcgcacgaggtcgggataacgcctgcaacgacatgggtaagatgctgcaattcgtgctgccgtggccacg
cagatccagcaggaggttatcaaaagcctatggcttcagctgcgacggggaagggttccttaagtgtgctcgttgcgaagtcctacgaagccaggaatcctgagatcgccagcctgtcaggcaagctgaaggcgtgtttctgccccatgaccctgccacccatggcctgctgctg
gtggcagcgtggccgctcctga