



Reference: PTXBC009938

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AAA-TOB3; TOB3; ATPase family AAA domain-containing protein 3B; AAA-ATPase TOB3; ATPase family, AAA domain containing 3B		

For research use only.

GENE DESCRIPTION	ATPase family, AAA domain containing 3B
GENE ID	83858
ORIGIN SPECIES	Human
ACCESSIONS	BC009938
VECTOR	pDONR223
NCBI SYMBOL	ATAD3B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]