

HLA-DPA1-major histocompatibility complex, class II, DP alpha 1 Gene

Reference: PTXBC009956

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MHC class II HLA-DPA1 antigen; DP(W3); DP(W4); HLA-DP1A; HLADP; HLASB; PLT1; HLA class II histocompatibility antigen, DP alpha 1 chain; HLA-SB alpha chain; MHC class II DP3-alpha; major histocompatibility complex, class II, DP alpha 1		

For research use only.

GENE DESCRIPTION	major histocompatibility complex, class II, DP alpha 1
GENE ID	3113
ORIGIN SPECIES	Human
ACCESSIONS	BC009956
VECTOR	pDONR223
NCBI SYMBOL	HLA-DPA1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcgccctgaagacagaatgtccatcatcagagctgtgatcttgagagccctctcttggcttctctgctgagtcctcaggagctggggccatcaaggcggaccatgtgtcaacttatgccgctttgtacagaccatagaccaacaggggagtttatgtttgaattgatgaagatgagca
gttctatgtggatctggataaaaaggagaccgtctggcatctggaggagtttggccgagcctttcttggagctcaggcgggctggctaaactgtctatattgaacaacaacttgaatacttgatccagcgttccaaccacactcaggccgccaatgatccccctgaggtgaccgtgttc
ccaaggagcctgtggagctgggagcccaacaccctcatctgccacattgacaggttcttccaccagtgtcaacgtcacatggctgtgcaatggggagccagtcactgaggggtgtcgtgagagccttcttcccagaacagattacagcttccacaagttccattacgtgaccttg
tgccctcagcagaggacgtctatgactgcaggggtggagcactgggcttggaccagccgtcctcaagcactgggagggccaagagccaatccagatgcctgagacaacggagactgtgctctgtgccctgggctgtgctggcctagtgggcatcatctgggaccgtcctcat
cataaagtctctgcgttctgccaatgacccccgggcccaggggcccctgtga