

SLC27A4-solute carrier family 27 (fatty acid transporter), member 4 Gene

Reference: PTXBC009959

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ACSVL4; FATP4; IPS; long-chain fatty acid transport protein 4; solute carrier family 27 (fatty acid transporter), member 4; solute carrier family 27 member 4		

For research use only.

GENE DESCRIPTION	solute carrier family 27 (fatty acid transporter), member 4
GENE ID	10999
ORIGIN SPECIES	Human
ACCESSIONS	BC009959
VECTOR	pDONR223
NCBI SYMBOL	SLC27A4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccctcagctgtctacgtctgcaaccggcgccgcatctggacggggcgccgcgccggagccgagccggcccaaatgctgttgagcctctctgtgggggtgctgtgttctccaagctggtgctgaaactgccctggaccaggtgggattctcctgttgttctcta
cttgggacttgccggtggcgcttcatccgggtcttcatgaagaccatcaggcctacttacttggtgatgtgctggtgatggagagctgggctacctgtactccgagaccgactggggacacgttccgctggaagggtgagaacgtgtccaccaccgaggtggaaggcacactcagc
cgctgtgacatggctgacgtggcgtgtatggtgctgaggtgccaggaaccgagggccggcggaatggctgctgtggccagcccaactggcaactgtgacctggagcgcttggctcaggtcttggagaaggaaactgccctgtatgcgcgcccatcttctgcgcctcctgcc
tgagctgcacaaacaggaaacctacaagttccagaagacagagctacggaaggaggccttgacccggcctattgtgaaagaccgctgttctatctagatgccagaaggccgctacgtcccgtggaccaagagcctacagccgcatccagcagcgaggagaagctgtga