

ARAP3-ArfGAP with RhoGAP domain, ankyrin repeat and PH domain 3 Gene

Reference: PTXBC009968

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CENTD3; DRAG1; arf-GAP with Rho-GAP domain, ANK repeat and PH domain-containing protein 3; ARF-GAP, RHO-GAP, ankyrin repeat and plekstrin homology domains-containing protein 3; Arf and Rho GAP adapter protein 3; PtdIns(3,4,5)P3-binding protein; centaurin-		

For research use only.

GENE DESCRIPTION	ArfGAP with RhoGAP domain, ankyrin repeat and PH domain 3
GENE ID	64411
ORIGIN SPECIES	Human
ACCESSIONS	BC009968
VECTOR	pDONR223
NCBI SYMBOL	ARAP3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaatgtggattggaccaccagcctctaaagccagcagcagtgaccagcagccagtggtctttacgacgccattctctctgacctggccgtcagaagtttggcactatgctttgctgctatccgtggggatgacagtgaggccaccctctcttgcacatcagaccctgcggc
gactacacaaccggaggacctgtccatgttcttccaatga