

ZNF581-zinc finger protein 581 Gene

Reference: PTXBC010097

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms:	HSPC189; zinc finger protein 581		

For research use only.

GENE DESCRIPTION	zinc finger protein 581
GENE ID	51545
ORIGIN SPECIES	Human
ACCESSIONS	BC010097
VECTOR	pDONR223
NCBI SYMBOL	ZNF581
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgtgtgtgtccatccccctgcctcagccttgccttctcctcgttgagacatggaggccctccccctcggacttgcctccccagaaactggacttctctctctccatcgatctccccaggttcatctctccaaggcccaaccactactgttattgacactcagggtgtccc
ctacacagtgtgtgtgacgaggtgacagagggagccaggccagtggtggtccaggccagaaaaagtgtctacagctgccccgtgtgtctcaagggtcttcgagtacatgtctacctcagcgacacagcatcaccactcggaggtaaagcccttcgagtgacatctgtggg
aaggcattcaagcgcgcagccacttggcacggcaccattcattcacttggcgggtgtgtggcgccccacggctgccccctctgccccctcgccttccgggatgcgggtgagctggccagcacagccgggtgcactctggggaacgcccgtttcagtgatccactgcctcgc
cgctttatggagcagaacactgcagaacacacgcggtggaagcatccatga