

C19orf10-chromosome 19 open reading frame 10 Gene

Reference: PTXBC010129

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0556 protein C19orf10; C19orf10; EUROIMAGE1875335; IL25; IL27; IL27w; R33729_1; myeloid-derived growth factor; interleukin 27 working designation; interleukin-25; stromal cell-derived growth factor SF20; myeloid derived growth factor		

For research use only.

GENE DESCRIPTION	chromosome 19 open reading frame 10
GENE ID	56005
ORIGIN SPECIES	Human
ACCESSIONS	BC010129
VECTOR	pDONR223
NCBI SYMBOL	C19orf10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcggccagcggagggtggaacggcgctggcgagcttggggccgctgctcctagggccgtggcgctgagggccggcgaggcgggtgcccagccccagacgggtggcggttgacgtgccccggcgggcgctgctgattccttccataacgtggggccgggg
gacaaatatacgtgtatgttcacttacgctctcaaggagggaaccaatgagcaatggcagatgagctggggaccagcgaagaccacagcacttcacctgcacatctggaggccccaggggaagtcctatctgtacttcacaggtcaaggcagaggtgcggggcgctgagattga
gtacgcatggcctactctaaagccgcatgtgaaagggaagtgatgccctctgaaaactgagggaattggaagtgaacaaacagcagtggtcagagggccggggcattcaaaagctgagctgtccaagctggtgattgtggccaaggcgcgcactgagctgtga