

FDX1L-ferredoxin 1-like Gene

Reference: PTXBC010155

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FDX2; ferredoxin-2, mitochondrial; adrenodoxin-like protein; adrenodoxin-like protein, mitochondrial; ferredoxin 1 like		

For research use only.

GENE DESCRIPTION	ferredoxin 1-like
GENE ID	112812
ORIGIN SPECIES	Human
ACCESSIONS	BC010155
VECTOR	pDONR223
NCBI SYMBOL	FDX1L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctagacatggccccctcctcaggagaactcgcggctgggctgccagattgtgctgacaccggagctggaaggagcgggaattcacctgcccaagatcaccaggaactctacgtggatgccatgtccccaagccccactga