

CLDN15-claudin 15 Gene

Reference: PTXBC010160

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	claudin-15; claudin 15		
For research use only.			

GENE DESCRIPTION	claudin 15
GENE ID	24146
ORIGIN SPECIES	Human
ACCESSIONS	BC010160
VECTOR	pDONR223
NCBI SYMBOL	CLDN15
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgtggtgtggaaccttggcttctcatggcaactgtgggtgctgtatgtgggtgactctgcaaacagctactggcgagtgctccactgtgcacgggaacgtcatcaccaccaacacatctcgagaacctctggttagctgtgccaccgactccctggcggtctacaact
gctgggagttccgtccatgctgtgacctctgtgctccagactccccagctcactgtcgggggactgtgtctctgtccggtgtgtctataaaggccctgtgaaggaggcgctcttgcaggttgagctccgtgtgaaggaggcggtgtgtgcaggggtatattc
aggcctgccgggcaactcatga