

PERP-PERP, TP53 apoptosis effector Gene

Reference: PTXBC010163

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PERP, TP53 apoptosis effector; KCPI; KRTCAP1; PIGPC1; THW; dJ496H19.1; p53 apoptosis effector related to PMP-22; 1110017A08Rik; KCP-1; keratinocyte-associated protein 1; keratinocytes associated protein 1; p53 apoptosis effector related to PMP22; p53-indu		

For research use only.

GENE DESCRIPTION	PERP, TP53 apoptosis effector
GENE ID	64065
ORIGIN SPECIES	Human
ACCESSIONS	BC010163
VECTOR	pDONR223
NCBI SYMBOL	PERP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatccgtgcggcctggcctgcgagcgtgccgtggatcctgccctgctcctactcagcgccatcgcttcgacatcagcgctggcggcggcggtggtgcagctagcgaccacggccagacgtcctcgtgtggtggaatgctcccaagaggcgggcagcgggt
cctacgaggagggtgtcagagcctcatggagtagcgtgggttagagcagcggtgccatgctctctgtggctcatcatcgtggtgatctgtttcatcctctctctcgcctctgtggacccagatgctgtctcctgagagtgtggaggtcctctgcttggctgtgttccag
atcatcctcctggttaattacccgtgaagtacaccagaccttcacccctcatgccaaccctgctgtcacttacatctataactgggcctacggcttgggtgggcagccagattatcgtgattggctgtgcttctctctgctgctcccaactacgaagatgaccttctgggcaatgccaa
gccaggtactctacacatctgctaa