

GINS2-GINS complex subunit 2 (Psf2 homolog) Gene

Reference: PTXBC010164

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSPC037; PSF2; Pfs2; DNA replication complex GINS protein PSF2; GINS complex subunit 2 (Psf2 homolog); GINS complex subunit 2		

For research use only.

GENE DESCRIPTION	GINS complex subunit 2 (Psf2 homolog)
GENE ID	51659
ORIGIN SPECIES	Human
ACCESSIONS	BC010164
VECTOR	pDONR223
NCBI SYMBOL	GINS2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacgctgccgaggtcgaaattcctgccgagaaaggagctggttaccattatcccaactcagctctggacaagatctacctcatcgggggggacgtgggcttttaacccctggtttaccctggaagtccccctgtggctggcgattaacctgaacaaagacagaaatgtcgcctgctccctccagagtgatgagtagagaaaagttggagaagatgagggatcatgaacgaaaggagaagaacttttaccctaatgccagcccttactacatggaacttacgaagctcctgttaaatcatgcttcagacaacatcccgaaggcagacgaaatccgacccctgtgcaaggatattgggacactcgtatagccaactccgagtgctgtgctgacagctttgtgagacagcaggaggacatgccaagctggataacttgacctgatggagatcaacaccagcgggactttctcacacaagcgtcaaccacatgtacaaactccgcgaacctccagcctctggagagtactcagctcaggactctag