

C14orf172-chromosome 14 open reading frame 172 Gene

Reference: PTXBC010167

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C14orf172; GCD14; Gcd14p; TRM61; hTRM61; tRNA (adenine(58)-N(1))-methyltransferase catalytic subunit TRMT61A; tRNA (adenine-N(1))-methyltransferase catalytic subunit TRM61; tRNA (adenine-N(1))-methyltransferase catalytic subunit TRMT61A; tRNA methyltran		

For research use only.

GENE DESCRIPTION	chromosome 14 open reading frame 172
GENE ID	115708
ORIGIN SPECIES	Human
ACCESSIONS	BC010167
VECTOR	pDONR223
NCBI SYMBOL	C14orf172
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagcttcgtggcatagcaggagctgatcaaggagggtgacacggccatctgtcactggccatgggtgcaatgggtggcagtgctgtgcagcgtggggcacagaccagaccggcatgggtgctctgcggcactcagttgacctatcgccgcccctcggtccaaggtgacgtg
cgcccgagggtggctgggtgtatgtgtctgcacccacgcccagctctggagcgtgaacctggcgaccgcagcagatcccttactccacagacatcgccctcatcaccatgatgttgagcttcggcccgctctgtgtgtgtgagctcggcaccggcagtggtctgtgtccacgc
catcatccgacattgcacccacgggtcactgcacacgtggaggtccaccagcagcgggcagagaaggccgggagggaggtccaggagcaccgtgtggccgctgggtgactgtgcgcaccaggacgtgtgccgagtggtttggcgtgagccacgtggccgacgcccgtct
tctggacatcccatcaccctggaggccgtggggccacgcctgggacgcctcaaggtcgaaggcggggcgttctgtctcttccacgtgcatcgagcaggtgcaacgcacatgccaggcgtggcagcgcgcggccttctcagagctgagcacctggaggtgctgccacaggtct
acaacgtgcgcactgtcagcctgccaccgcccagctgggacaggcacagatggccctgcccgtctcgacaccagccccttcgcagcggcagcgcctatgaaggaggccgtggggccacaccggctacctgaccttcggccacaagaccagcagctag