

C9orf85-chromosome 9 open reading frame 85 Gene

Reference: PTXBC010179

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C9orf85; chromosome 9 open reading frame 85		

For research use only.

GENE DESCRIPTION	chromosome 9 open reading frame 85
GENE ID	138241
ORIGIN SPECIES	Human
ACCESSIONS	BC010179
VECTOR	pDONR223
NCBI SYMBOL	C9orf85
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagctcccagaaggcaacgtggctcgttcagacctcagaagcaccagaatacgttttagcttcaaaatgacaagttcgataaaagtgtcagaccaagtgttaa