

DPH3-DPH3, KTI11 homolog (S. cerevisiae) Gene

Reference: PTXBC010181

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	DPH3, KTI11 homolog (S. cerevisiae)
GENE ID	285381
ORIGIN SPECIES	Human
ACCESSIONS	BC010181
VECTOR	pDONR223
NCBI SYMBOL	DPH3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagtggtttcatgacgaggtggaaatcgaggactccaatgatgacgaggactcggagacgtatttctatccctgccatgtggagataacttctccatcaccaaggagatttggagaatggggaagacgtggcaacgtgtcctagctgctctctcattataaaagtattatgacaaaga
tcagtttgtgtggagaaacagtcgccagccccctcagccaacaaagaattagttaaatgctga