

RPS4Y1-ribosomal protein S4, Y-linked 1 Gene

Reference: PTXBC010286

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RPS4Y; 40S ribosomal protein S4, Y isoform 1; 40S ribosomal protein S4, Y; ribosomal protein S4Y; ribosomal protein S4, Y-linked 1		

For research use only.

GENE DESCRIPTION	ribosomal protein S4, Y-linked 1
GENE ID	6192
ORIGIN SPECIES	Human
ACCESSIONS	BC010286
VECTOR	pDONR223
NCBI SYMBOL	RPS4Y1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccggggccccaagaagcacttaagcgtgtgtcagcgccgaagcattggatgcttgacaaactaacgggtgtattgacatcgtccatcgacagggtccccaagaagctgagggaatgtcttctctgacgtcttctcaggaatagactcaagtagcgttgactggagatgaggta
aagaagatatgtatgaacgtttcatcaaaattgatggcaaggttcgagtgatgtcacatccctgctggattcattggtatgcacagcagcagagaagacaggtgaacatttccgctgtctatgacccaaggccggtttgtctgttaccgcatcacagtggaaaggcaagtagacaagt
gtgcaagtgagggaagattactgtgggagtgaaagggaatccctcacctggtgactcatgatgctcgaaccatccgctaccagatcctgtcatcaaggtgaacgatactgtgcagattgatttagggactggcaagataatcaactttatcaaaattgatacaggcaatttgtgtatggtgattgg
tggagccaacctcgtcgtgtgtgtgtatcaccaacagggaagacatcctggtctttgatgtggtgcatgtgaaggatgccaatggcaacagcttggccagaggctttcaacattttgtcattggcaatggcaataaaccttggatttcctgccaggggaaaggcattcgacttac
tgttctgaagagagagataagagcctggccaccaaacagagcagtggtgctaa