

RNF111-ring finger protein 111 Gene

Reference: PTXBC010369

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ARK; E3 ubiquitin-protein ligase Arkadia; ring finger protein 111		
For research use only.			

GENE DESCRIPTION	ring finger protein 111
GENE ID	54778
ORIGIN SPECIES	Human
ACCESSIONS	BC010369
VECTOR	pDONR223
NCBI SYMBOL	RNF111
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggttctgatatggcaggctatcctcacatcgttacatttcacaggatggatggaacatcagaggtccttcaggggcaattttgaggaaactgattcatttggagaaagattaggcaatgtcaatcgtggagcatccaggggacaattgaaagatgtacatatccataaataca
aaaagaggaaactgcactgcaacaagatggggaagggaaggactgaggaaacacagaggaaaaatgtactatctgtttgtctattttagaggaaagtgagatgtgagacgtcttccatgtatgcacctttccaccaagtgtgtgtgaccaatggttgattaccaataagaagtgcctcat
atgcagagtgacattgaggccagctgccaagtgaagtga