

CSTA-cystatin A (stefin A) Gene

Reference: PTXBC010379

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AREI; PSS4; STF1; STFA; cystatin-A; cystatin A (stefin A); cystatin AS; cystatin A		

For research use only.

GENE DESCRIPTION	cystatin A (stefin A)
GENE ID	1475
ORIGIN SPECIES	Human
ACCESSIONS	BC010379
VECTOR	pDONR223
NCBI SYMBOL	CSTA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatacctggaggcttatctgaggccaaacccgactccagaatccaggagattgttgataaggttaaccacagcttgaagaaaaacaaatgagacttatgaaaattggaagctgtgcagtataaaactcaagttgttgcctggaacaaattactacattaggtacgagcaggtgata
ataaatatgacacttgaagattcaaaagtctcccgacaaaatgaggacttggacttactggataccaggttgacaaaaacaggatgacgagctgacggccttttag