



Reference: PTXBC010381

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HPR1; P84; P84N5; THO complex subunit 1; hTREX84; nuclear matrix protein p84; tho1; THO complex 1		

For research use only.

GENE DESCRIPTION	THO complex 1
GENE ID	9984
ORIGIN SPECIES	Human
ACCESSIONS	BC010381
VECTOR	pDONR223
NCBI SYMBOL	THOC1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]