

GNG10-guanine nucleotide binding protein (G protein), gamma 10 Gene

Reference: PTXBC010384

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	guanine nucleotide-binding protein G(I)/G(S)/G(O) subunit gamma-10; guanine nucleotide binding protein (G protein), gamma 10; guanine nucleotide binding protein 10; G protein subunit gamma 10		

For research use only.

GENE DESCRIPTION	guanine nucleotide binding protein (G protein), gamma 10
GENE ID	2790
ORIGIN SPECIES	Human
ACCESSIONS	BC010384
VECTOR	pDONR223
NCBI SYMBOL	GNG10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctccggggctagcgcgagcgcctgcagcgttgtagagcagctcaagttggagcgtggcgtggagagatcaaggtctctcaggcagctgcagagcttcaacagtactgtatgcagaatgcctgcaaggatccctgctgtgggtgttccagctggaagtaacccctccg
ggagcctagatcctgtgctttactctga