

SOCS2-suppressor of cytokine signaling 2 Gene

Reference: PTXBC010399

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CIS2; Cish2; SOCS-2; SSI-2; SSI2; STAT2; suppressor of cytokine signaling 2; CIS-2; STAT-induced STAT inhibitor 2; cytokine-inducible SH2 protein 2		

For research use only.

GENE DESCRIPTION	suppressor of cytokine signaling 2
GENE ID	8835
ORIGIN SPECIES	Human
ACCESSIONS	BC010399
VECTOR	pDONR223
NCBI SYMBOL	SOCS2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaccctgcgggtgccttgagccctccgggaatggcggggaaggacgggagccagtgggggaccgcggtcgccggagagccatcccgaggcggcggtctggcgaaggccctgcgggagctcggtcagacaggatggtagtgggaagtatgactgttaatgaagc
caaagagaaataaaaggccaccagaaggaaacttcttgattagagatagctcgattcagactactactaacaatatctgttaaaacatcagctggaccaactaatcttcgaatcgaataccaagacggaaaattcagattggactctatcatatgtgtcaaatccaagcttaacaatttgac
agtggtgtcatctgactactatgttcagatgtgcaaggataagcggacaggtccagaagccccgggaacgacgtgttcacctttatctgacaaaaccgtctacacgtcagcaccatctctgcagcatctctgtaggtcaccattacaatgtaccgggtccatctggggactgc
ctttaccacaagactaaaagattcttggaagaatataaattccaggtataa