



Reference: PTXBC010423

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PVRL4; EDSS1; LNIR; PRR4; nectin-4; Ig superfamily receptor LNIR; poliovirus receptor-related 4; poliovirus receptor-related protein 4; nectin cell adhesion molecule 4		

For research use only.

GENE DESCRIPTION	poliovirus receptor-related 4
GENE ID	81607
ORIGIN SPECIES	Human
ACCESSIONS	BC010423
VECTOR	pDONR223
NCBI SYMBOL	PVRL4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]