

C14orf179-chromosome 14 open reading frame 179 Gene

Reference: PTXBC010436

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C14orf179; CED3; intraflagellar transport protein 43 homolog; IFT complex A subunit; intraflagellar transport 43 homolog; intraflagellar transport 43		

For research use only.

GENE DESCRIPTION	chromosome 14 open reading frame 179
GENE ID	112752
ORIGIN SPECIES	Human
ACCESSIONS	BC010436
VECTOR	pDONR223
NCBI SYMBOL	C14orf179
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggatttgctcacttggacaggagcttcctacagcttgctacctccaggcccaagatgggtgcggagctcaacaggagtcagcgcaggccgagaatcacctcaatggcaagaattcctttgactctgactggagagacttctctgctaaattacctcgtccgacagg
gaggctgggcagggtgattccgtgaaggcttcgaacgggtacccaaacaggcaacaacagctggatctgaacgcatgctatcacaaaacgcatcacagaaattgggctggttcattggaaaggcagatattcctatcattccggatctggaggagtagaggaagactttgtttg
cagggtggcagccccccagcatccagataaagcgggtgatgacctaccgtgacctggacaatgacctcatgaagtactcagccattcagacactggatggggagatcgacctgaaactcctcacaaagtctcgcgcgggagcacgaagtcgggaggtgatgtcggctgggact
gggaccatctgttactgaggtgtcctcagaggctcctcactgagtgggacccactgcagacggagaaggagaccctgcggggcaggccaggcacacctga